


Pilon Fusion System

Surgical Technique



Introduction

The Arthrex Pilon Fusion System is designed to treat distal tibia fractures that require both fracture reduction and primary ankle arthrodesis. Severe damage to the tibiotalar (TT) joint often results in posttraumatic arthritis, pain, stiffness, and the need for secondary surgeries. The Primary Pilon Fusion System provides another option to address these severe fracture patterns with primary TT arthrodesis of the articular surface to help avoid secondary surgery and chronic pain. Anterolateral and posterior approaches, depending on the fracture pattern, allow for fracture management, lengthy bridging techniques, anatomic implants, and fracture-specific locking configurations.

Arthrex offers comprehensive solutions to treat these patients with the ArthroFX® large external fixation system, Ankle Fusion Plating System, Titanium Ankle Fracture Management System, Distal Tibia Plating System, and the FibuLock® fibular nail.

Note: The Ankle Fusion System is required for Pilon Fusion system use




Anterolateral Pilon/Tibiotalar Fusion Plate Design

The tibial shaft and talus fixation points can use a combination of 4.5 mm locking, 4.5 mm cortical, and 5.5 mm cancellous screws for increased strength across the ankle joint. The anterolateral plate has two rows of distal 3.5 mm locking or 3.5 mm cortical screws to address complex, high-energy pilon fractures. The use of 3.5 mm screws distally allows for a high density of screw fixation and additional options for fracture reduction. An oblong slot and K-wire/BB-Tak holes facilitate proper plate placement and provisional fixation. An oblong slot and K-wire/BB-Tak holes facilitate proper plate placement and provisional fixation.

Note: 15-hole and 18-hole plates are only available in sterile.

The anatomic distal contour allows for robust fixation across the ankle joint while providing the maximum number of fixation points for a successful fusion.


- Plate lengths: 5-hole (112 mm), 7-hole (137 mm), 9-hole (163 mm), and 12-hole (201 mm)
- Sterile options: 15-hole (239 mm) and 18-hole (277 mm)




Anterolateral Primary Pilon Fusion Plate Technique



Use an anterior or anterolateral approach, ensuring the neurovascular structures are protected. The fracture can be addressed with provisional fixation to aid in restoring the articular surface of the distal tibia.



Fracture management and joint preparation are managed accordingly based on the fracture pattern and surgeon preference. The ankle joint can be prepped for a fusion by using the various chisels and curettes offered in the system.




Attach the 4.5 mm percutaneous insertion handle to the most distal shaft screw of the plate. Slide the plate proximally in the previously created submuscular pocket. It is recommended to slide the plate proximally first and then seat it distally into the appropriate position based on anatomic contours.

The percutaneous insertion handle is cannulated to enable provisional fixation with a K-wire if desired.


Anterolateral Primary Pilon Fusion Plate Technique (Cont.)



Once the plate is positioned, drill with the 3.0 mm drill bit through the 3.0 mm/4.5 mm drill guide. Measure the appropriate screw length with the depth gauge and place a 4.5 mm cortical screw in the oblong slot; conversely, a BB-Tak can be used to provide provisional fixation.




Provisionally fix the talus with a BB-Tak. Drill with the 3.0 mm calibrated drill through the 4.5 mm locking drill guide; the screw length can be measured off the back of the drill guide. Using the T20 driver, insert the appropriate locking or cortical screw.



Use 1.6 mm K-wires to secure the plate to the distal tibia and to secure the fracture fragments. Use a 3.5 mm cortical screw to manipulate the bone toward the plate or a 3.5 mm locking screw to maintain spatial positioning.


Drill with the 2.5 mm calibrated drill through the 3.5 mm locking drill guide. Measure with the drill guide or the depth gauge and implant a 3.5 mm cortical or locking screw with the T15 driver.

Anterolateral Primary Pilon Fusion Plate Technique (Cont.)



Once the distal fixation is completed, add a 4.5 mm locking screw in the kickstand screw hole to provide additional stability across the joint. Use the 4.5 mm locking drill guide drill with the 3.0 mm calibrated drill to measure the screw length off the drill guide.


Note: A 4.5 mm nonlocking cortical or 5.5 mm cancellous screw may be used if desired.



Distal screw trajectories.



Once the distal fixation is completed, add 4.5 mm cortical nonlocking or locking screws proximally as needed.



Optional: AlloSync™ demineralized bone graft hydrated with concentrated platelet-rich plasma (PRP) from bone marrow aspirate (BMA) can be used to augment the ankle arthrodesis procedure.


Posterior Pilon/Tibiotalar Fusion Plate Design

The posterior plate has two rows of distal 4.5 mm locking or 4.5 mm nonlocking cortical screws to address complex, high-energy pilon fractures while providing stability and strength for a successful fusion. With three points of fixation in the talus, along with the locking kickstand screw across the ankle joint, the posterior plate provides a robust buttress effect for addressing fracture reduction while providing the spatial arrangement and stability needed to maintain talus positioning for tibiotalar arthrodesis.


An oblong slot and K-wire/BB-Tak holes facilitate proper plate placement and provisional fixation. The anatomic distal contour allows for robust fixation across the ankle joint while providing the maximal number of fixation points.

- Plate lengths: 5-hole (119 mm), 7-hole (144 mm), 9-hole (170 mm), 12-hole (208 mm)
- Sterile options: 15-hole (246 mm)

Note: 15-hole plate is only available in sterile.



Accepts 4.5 mm cortical, 4.5 mm locking, and 5.5 mm cancellous screws



Posterior Pilon/Tibiotalar Fusion Plate Technique




Use a posterolateral or posteromedial approach, dictated by fracture pattern and surgeon preference, while carefully protecting tendons and neurovascular structures. Obtain provisional fracture reduction and prepare the ankle joint surfaces for fusion.




Attach the 4.5 mm percutaneous insertion handle to the most distal shaft screw of the plate. Slide the plate proximally past the ideal insertion point and then move it distally until it seats in the appropriate position on the distal tibia poster facet of the talus.

The percutaneous insertion handle is cannulated to enable provisional fixation with a K-wire if desired.


Posterior Pilon/Tibiotalar Fusion Plate Technique (Cont.)




Place a BB-Tak in the oblong hole for provisional fixation and begin fixation with 4.5 mm cortical screws. Drill with the 3.0 mm drill bit through the 3.0/4.5 mm drill guide.



Measure the length with the depth gauge and insert the desired 4.5 mm screw. Remove the BB-Tak.




Add additional 4.5 mm screws as needed. Strategically placed screws may help with manipulations of the fracture fragment position.




Secure the talus by placing a BB-Tak in the medial or lateral talar screw hole. A 4.5 mm cortical screw can help decrease any gaps between the bone and the plate, but may affect alignment. Be sure to maintain the desired foot and fracture alignment while drilling. The use of 4.5 mm locking screws will help maintain fracture alignment. Drill with the 3.0 mm drill bit and the 3.0 mm/4.5 mm drill guide.

Posterior Pilon/Tibiotalar Fusion Plate Technique (Cont.)




Once the distal fixation is completed, add a 4.5 mm locking screw to the kickstand screw hole to provide additional stability across the joint. Use the 4.5 mm locking drill guide drill with the 3.0 mm calibrated drill to measure the screw length off the drill guide.

Note: A 4.5 mm nonlocking cortical or 5.5 mm cancellous screw can also be used.




Distal screw trajectories.

Note: The locking screw trajectories of the kickstand screw and the central talar screw will intersect. The surgeon must be mindful during drilling and select screws of the appropriate length to avoid interference.




Once the distal fixation is completed, add a 4.5 mm cortical or locking screw as needed.




Optional: AlloSync™ demineralized bone graft hydrated with concentrated PRP from BMA can be used to augment the ankle arthrodesis procedure.


Pilon Fusion Tray Overview



Level 1
Implant Tray



Level 2
Instrument Tray



Screw Caddy
3.5 mm Locking and
Cortical Screws

Note: Ankle Fusion System is required for Pilon Fusion System use.

Supporting Products


Arthrex Ankle Fusion Plating System

The titanium Ankle Fusion Plating System provides a complete solution for ankle fusion management, with a comprehensive offering of anatomy-specific plates available for either tibiotalar or tibiotalocalcaneal arthrodesis. A variety of screw options—including locking, nonlocking, cortical, cancellous, and hybrid designs—are provided to address all fixation needs. Specific instrumentation designed to help gain access to and prepare the fusion sites is also included in the set. The Ankle Fusion Plating System was designed to provide the solution to your ankle fusion fixation needs.


- Comprehensive instrumentation for joint preparation, distraction, and compression and assistance with optimal fixation
- Maximum fixation points within each plate

System Features


- Anatomically designed for use with three surgical approaches: anterior, lateral, posterior
- Four compression modes available in system
 - Anatomic compression hole
 - Oblong compression hole
 - Mini joint compressor/distractor
 - 6.7 mm cannulated lag screws or 7.0 mm XL Compression FT screws




Note: Ankle Fusion System is required for Pilon Fusion System use.



Anterior Tibiotalar Plate
AR-8970AR



Lateral Tibiotalar Plate
AR-8970TT



Lateral Tibiotalocalcaneal
Plate – AR-8970TTC



Posterior Tibiotalocalcaneal
Plate – AR-8970PR

Supporting Products – Ankle Fusion Plating System

Cannulated Screws

The versatility of the Ankle Fusion System provides a complete solution for treating ankle arthritis in one comprehensive instrument case. The instrument set can be configured to house either 6.7 mm cannulated lag screws or 7.0 mm XL Compression FT screws for percutaneous compression across the arthrodesis site.




6.7 mm LPS Screw (18 mm thread)



7.0 mm XL Compression FT Screw

- **Low-profile head** – 1.0 mm shorter than a traditional 6.5 mm AO screw, while still using a 3.5 mm hex driver.
- **Increased pull-out** – 30% better than a standard 6.5 mm AO screw.¹
- **Deeper threads** – Using a 2.4 mm guide pin allows the threads to be deeper than a standard AO screw.
- **Self-drilling/tapping** – Speeds up the insertion process.
- **Headless design** – Minimal risk of impingement or soft-tissue irritation.
- **Fully threaded compression** – Variable-stepped thread pitch and tapered proximal profile work together to compress bone fragments with the purchase of a fully threaded screw.
- **Self-drilling/tapping** – Helical relief flutes assist in bone removal to reduce insertion torque.



Osteotome, angled up, 5.5 mm

Osteotome, angled up, 12 mm

Osteotome, straight, 5.5 mm

Osteotome, straight, 12 mm

Cup Curette, straight, 6.0 mm

Cup Curette, curved, 6.0 mm

Joint Preparation


Straight and curved curettes and osteotomes have been added to the ankle fusion tray to help with the removal of cartilage from the ankle and subtalar joints. These instruments are appropriately designed for the ankle and come standard in each ankle fusion tray, simplifying joint preparation in the OR setting.

Mini Joint Compressor/Distractor


Adaptable for distraction and compression of arthrodesis sites, this unique device facilitates joint preparation and allows for excellent compression prior to definitive fixation. The device uses 1.6 mm or 2.4 mm guidewires or 3.0 mm traction screws, which are included in the system.




Ankle Fusion Tray Overview




Level 1
Instrument Case
Implant Tray and
Instruments
AR-8970C-01




Level 2
Interchangeable
Cannulated Screw
Instrument Tray




Level 2
7.0 mm Screw Instrumentation
■ AR-8970C-70 – 7.0 mm Tray
■ AR-8970C-SC-70 – 7.0 mm
Screw Caddy Insert




Level 2
6.7 mm Screw Instrumentation
■ AR-8970C-67 – 6.7 mm Tray
■ AR-8970C-SC-67 – 6.7 mm
Screw Caddy Insert



Level 3
 Auxiliary
 Instruments and
 Screw Caddy



Screw Caddy
 Interchangeable
 Cannulated Screws



Screw Caddy Inserts

- 7.0 mm
 Compression FT
 Screw Caddy
 (AR-8770-SC-70)
- 6.7 mm Cannulated
 Screw Caddy
 (AR-8970C-SC-67)

Supporting Products

Biologic Options

Angel® Concentrated PRP System


Technology is what sets the Angel system apart from the competition. The Angel system is the only one to provide PRP concentrate from BMA with adjustable cellular levels. Bone marrow is a rich source of platelets and nucleated and progenitor cells. Customization of cellular levels is necessary to reduce the number of neutrophils in BMA, which can be detrimental to bone healing.

Features and Benefits


- Proprietary platelet sensor system
- Adjustable platelet concentrations
- Adjustable WBC concentrations
- Flexible processing volume 40 mL to 180 mL
- Each processing kit can process 3 cycles up to 180 mL on the same patient
- Programmable – can store up to 30 custom processing protocols
- Closed system: delivers PRP, platelet-poor plasma, and RBCs into separate, sterile compartments



In vitro culture expansion of progenitor cells²



48 hours



96 hours

| Angel cPRP System | Platelet Concentration (K/ μ L) | Nucleated Cell Concentration (K/ μ L) | Hematopoietic Cell Concentration (K/ μ L) | Total Neutrophils ($\times 10^6$) |
|--------------------------|-------------------------------------|---|---|-------------------------------------|
| BMA | 87.7 \pm 6.4 | 24.5 \pm 15.6 | 0.002 \pm 0.001 | 612.1 |
| PRP Concentrate From BMA | 787.0 \pm 317.6 | 240.5 \pm 186.6 | 0.081 \pm 0.056 | 132.9 |
| Increase Above Baseline | $\sim 9\times$ | $\sim 10\times$ | $\sim 33\times$ | 80% \downarrow |

Data from Arthrex, Inc. Data on file (APT-02569). Naples, FL; 2018.

AlloSync™ Pure Demineralized Bone Matrix


AlloSync Pure demineralized bone matrix (DBM) is derived from 100% human allograft bone with no extrinsic carriers. AlloSync Pure bone matrix resists irrigation and can be used in a fluid environment. The clinician can control the handling properties of AlloSync Pure bone matrix, which includes decreasing the viscosity for injectable applications or increasing the viscosity for open procedures. The proprietary rice-shape fiber technology used to process AlloSync Pure bone matrix increases the osteoinduction and osteoconductive surface area to accelerate cellular ingrowth.³



Supporting Products – Biologic Options

AlloSync™ Demineralized Bone

Demineralized cancellous sponges and cortical fibers are optimal for combination with blood products such as concentrated BMA. When combined with BMA, AlloSync demineralized bone grafts provide an optimal environment for bone formation that includes the “healing triad”: cell, signal and scaffold.⁴



Angel® System

| Product Description | Item Number |
|----------------------------------|-------------|
| Angel System | ABS-10060 |
| Angel Bone Marrow Processing Kit | ABS-10062 |
| Angel Blood Access Kit | ABS-10067 |

AlloSync Cancellous Sponges


| Product Description | Item Number |
|--------------------------------|-------------|
| Cube, 8.0 mm × 8.0 mm × 8.0 mm | ABS-2005-01 |
| Cube, 10 mm × 10 mm × 10 mm | ABS-2005-02 |
| Cube, 12 mm × 12 mm × 12 mm | ABS-2005-03 |
| Strip, 10 mm × 10 mm × 3.0 mm | ABS-2006-01 |
| Strip, 15 mm × 40 mm × 3.0 mm | ABS-2006-02 |
| Strip, 26 mm × 19 mm × 7.0 mm | ABS-2006-03 |
| Strip, 10 mm × 20 mm × 7.0 mm | ABS-2006-04 |
| Chips (1.0 mm-4.0 mm), 1.0 cc | ABS-2007-01 |
| Chips (1.0 mm-4.0 mm), 2.5 cc | ABS-2007-02 |
| Chips (1.0 mm-4.0 mm), 5.0 cc | ABS-2007-03 |

AlloSync Cortical Fibers

| Product Description | Item Number |
|---------------------|-------------|
| Fibers, 1.0 cc | ABS-2008-01 |
| Fibers, 2.5 cc | ABS-2008-02 |
| Fibers, 5.0 cc | ABS-2008-03 |
| Fibers, 10 cc | ABS-2008-04 |


AlloSync Pure DBM

| Product Description | Item Number |
|---------------------|-------------|
| Pure DBM, 1.0 cc | ABS-2010-01 |
| Pure DBM, 2.5 cc | ABS-2010-02 |
| Pure DBM, 5.0 cc | ABS-2010-05 |
| Pure DBM, 10 cc | ABS-2010-10 |



Supporting Products

Snap-Off Compression FT Pins



The new Snap-Off Compression FT pins offer solutions for upper and lower extremities and the trauma setting. The ease of use, variable-stepped thread pitch, and wide range of sizes in diameters of both 1.9 mm and 2.4 mm represent an innovative new technology to help treat pathologies including joint arthrodesis, hammertoe, intra- and extra-articular fractures, and nonunions. The system can be used in conjunction with other Arthrex products to compress bone fragments.

- **Variable-stepped thread pitch**—the screw tip’s wider thread pitch enters the bone faster than trailing threads, gradually compressing the fragments as the screw is advanced
- Convenient and quick surgical technique
- Self-drilling and self-tapping to facilitate insertion
- Multiple size options available in both 1.9 mm and 2.4 mm diameters
- Available in titanium

Snap-Off Compression FT Pins (AR-9938S)


| Product Description | Item Number |
|--|---------------|
| Instruments | |
| Pin Insertion Guide | AR-9938-01 |
| Measuring Drill Guide | AR-9938-02 |
| Depth Device | AR-9938-03 |
| Locking Pin Driver | AR-9938-09 |
| Pin Cutter | AR-9938-11 |
| Cannulated Scalpel Handle, 1.35 mm | AR-9938-12 |
| Implants | |
| 1.9 mm Snap-Off Compression FT Pins, 10 mm-50 mm (2.0 mm increments) | AR-9919-10-50 |
| 2.4 mm Snap-Off Compression FT Pins, 10 mm-50 mm (2.0 mm increments) | AR-9924-10-50 |
| Compression FT Pin, solid (sterile), 1.9 mm × 30 mm | AR-9919T-30S |
| Compression FT Pin, solid (sterile), 1.9 mm × 50 mm | AR-9919T-50S |
| Compression FT Pin, solid (sterile), 2.4 mm × 30 mm | AR-9924T-30S |
| Compression FT Pin, solid (sterile), 2.4 mm × 50 mm | AR-9924T-50S |
| Disposables for 1.9 mm Snap-Off Compression FT Pins | |
| Calibrated Measuring K-Wire, 1.1 mm, qty. 6 | AR-9938-04 |
| Hard-Bone Drill, 1.5 mm, qty. 2 | AR-9938-05 |
| Pin Extractor, 1.9 mm, qty. 2 | AR-9938-08 |
| Disposables for 2.4 mm Snap-Off Compression FT Pins | |
| Calibrated Measuring K-Wire, 1.35 mm, qty. 6 | AR-9938-10 |
| Hard-Bone Drill, 1.7 mm, qty. 2 | AR-9938-06 |
| Pin Extractor, 2.4 mm, qty. 2 | AR-9938-07 |

Ordering Information

Primary Pilon Fusion System (AR-8963FS)

| Product Description | Item Number |
|--|---------------|
| Screwdriver, T15 hexalobe | AR-8943-10 |
| Drill Guide, 3.5/2.5 mm | AR-8943-14 |
| Drill Guide, locking, 3.5 mm | AR-8943-43 |
| Countersink, 3.5/4.0 mm | AR-8950-03 |
| Depth Device, 3.5/4.0 mm | AR-8963-13 |
| Chisel, straight, 14 mm | AR-8963-14 |
| Chisel, curved, 14 mm | AR-8963-15 |
| Chisel, curved, 7.0 mm | AR-8963-16 |
| Chisel, straight, 7.0 mm | AR-8963-17 |
| Cup Curette, 9.0 mm | AR-8963-18 |
| Lag Drill Guide, 3.5 mm/2.5 mm | AR-8963-21 |
| Percutaneous Insertion Handle, 4.5 mm | AR-8963-24 |
| Driver, T15 hexalobe, 6.0 in, AO, qty. 2 | AR-8963-25 |
| Driver, T15 hexalobe, 6.0 in, straight, AO, qty. 2 | AR-8963-26 |
| Primary Pilon Fusion Case | AR-8963C-03 |
| Plates | |
| Anterolateral Pilon Fusion Plate, 5h, left | AR-8963ALF-05 |
| Anterolateral Pilon Fusion Plate, 7h, left | AR-8963ALF-07 |
| Anterolateral Pilon Fusion Plate, 9h, left | AR-8963ALF-09 |
| Anterolateral Pilon Fusion Plate, 12h, left | AR-8963ALF-12 |
| Anterolateral Pilon Fusion Plate, 5h, right | AR-8963ARF-05 |
| Anterolateral Pilon Fusion Plate, 7h, right | AR-8963ARF-07 |
| Anterolateral Pilon Fusion Plate, 9h, right | AR-8963ARF-09 |
| Anterolateral Pilon Fusion Plate, 12h, right | AR-8963ARF-12 |
| Posterior Pilon Fusion Plate, 5h | AR-8963PF-05 |
| Posterior Pilon Fusion Plate, 7h | AR-8963PF-07 |
| Posterior Pilon Fusion Plate, 9h | AR-8963PF-09 |
| Posterior Pilon Fusion Plate, 12h | AR-8963PF-12 |

| Product Description | Item Number |
|--|----------------|
| Sterile Plates | |
| Anterolateral Pilon Fusion Plate, 15h, left | AR-8963ALF-15S |
| Anterolateral Pilon Fusion Plate, 18h, left | AR-8963ALF-18S |
| Anterolateral Pilon Fusion Plate, 15h, right | AR-8963ARF-15S |
| Anterolateral Pilon Fusion Plate, 18h, right | AR-8963ARF-18S |
| Posterior Pilon Fusion Plate, 15h | AR-8963PF-15S |
| 3.5 mm Screws, Low Profile, Ti | |
| 3.5 mm × 20 mm-60 mm (2.0 mm increments) | AR-8935-20-60 |
| 3.5 mm × 65 mm-80 mm (5.0 mm increments) | AR-8935-65-80 |
| 3.5 mm × 90 mm-120 mm (10 mm increments) | AR-8935-90-120 |
| 3.5 mm Screws, Low Profile, Ti, locking | |
| 3.5 mm × 20 mm-50 mm (2.0 mm increments) | AR-8935L-20-50 |
| 3.5 mm × 55 mm-60 mm (2.0 mm increments) | AR-8935L-55-60 |
| Disposables | |
| Drill Bit, 3.5 mm, qty. 2 | AR-4160-35 |
| Drill Bit, 2.5 mm, qty. 2 | AR-8963-19 |
| Guidewire, drill tip, 3.0 mm, qty. 4 | AR-8963-20 |
| Drill Bit, calibrated, long, 3.0 mm | AR-8970-30L |
| Bone Tap | AR-8963-23 |



Ankle Fusion Plating System, 7.0 mm Set (AR-8970S-70)

| Product Description | Item Number |
|---|----------------|
| Perc Drill Guide, Compression FT | AR-8750-02 |
| Drill Guide, threaded, locking, 4.5 mm, qty. 2 | AR-8970-01 |
| Drill Guide, 3.0 mm/4.5 mm | AR-8970-02 |
| Drill Guide, 3.0 mm/5.5 mm | AR-8970-05 |
| Depth Measuring Device, long, 4.5 mm/5.5 mm | AR-8970-07L |
| Depth Device, cannulated screws | AR-8750-01 |
| Drive Shaft, T20 hexalobe, qty. 2 | AR-8970-03 |
| Driver, T20 hexalobe, straight | AR-8970-04 |
| Driver, T20 hexalobe, straight, AO, qty.2 | AR-8970-08 |
| Driver, T25 hexalobe, ISO, cannulated, qty. 2 | AR-8770-01 |
| Driver, T25 hexalobe, ISO, solid | AR-8770-04 |
| Ratcheting Handle, cannulated, large AO handle, QC | AR-8970RH |
| Mini Joint Distractor/Compressor | AR-8970JD |
| Axial Handle, trilobe QC, ratcheting | AR-8770RH |
| Soft Tissue Protector, 2.4 mm | AR-8770-06 |
| Bone Reduction Forceps, qty. 2 | AR-8943-07 |
| Hohmann Retractor, 9.5 in, 17 mm pointed, qty. 2 | AR-9260-34 |
| Cup Curette, straight shaft, 6.0 mm | AR-8970-11 |
| Cup Curette, curved shaft, 6.0 mm | AR-8970-12 |
| Cobb Elevator, 9.0 mm | AR-8640 |
| Small Joint Osteotome Angled Up, 0.217 in (5.5 mm), w/ handle | AR-8650-08 |
| Small Joint Osteotome Straight, 0.217 in (5.5 mm), w/ handle | AR-8650-09 |
| Small Joint Osteotome Angled Up, 0.472 in (12 mm), w/ handle | AR-8970-13 |
| Small Joint Osteotome Straight, 0.472 in (12 mm), w/ handle | AR-8970-14 |
| Screw Holding Forceps | AR-8941F |
| Guidewire Sleeve Insert, 1.6 mm | AR-8970-06 |
| Ankle Fusion Plating System Instrument Case | AR-8970C-01 |
| Ankle Fusion Instrument Case, 7.0 mm tray | AR-8970C-70 |
| Ankle Fusion Caddy, 7.0 mm insert | AR-8970C-SC-70 |

Disposables for AR-8970S-70 (not included in set, order separately)

| Product Description | Item Number |
|-------------------------------------|--------------|
| Drill Bit, calibrated, long, 3.0 mm | AR-8970-30L |
| Drill Bit, cannulated, long, 3.0 mm | AR-8970-30CL |
| Drill Bit, cannulated, long, 4.5 mm | AR-8970-45CL |
| Drill Bit, long, 4.5 mm | AR-8970-45L |
| Drill Bit, long, 5.5 mm | AR-8970-55L |
| Drill Bit, cannulated, long, 5.5 mm | AR-8970-55CL |
| Drill Bit, cannulated, 5.0 mm | AR-8770-02 |
| Profile Drill, X-large, 7.0 mm | AR-8770-03 |
| BB-Tak, large | AR-8970-09 |
| BB-Tak, large, threaded | AR-8970-09T |
| Traction Screw, 20 mm | AR-8950JD-2 |

Disposables for AR-8970S-70 (not included in set, order separately)

| Product Description | Item Number |
|--|---------------|
| Traction Post, threaded, 4.5 mm | AR-8970JD-45S |
| Guidewire w/ Trocar Tip, 0.095 in (2.4 mm) × 9.25 in | AR-8770K |
| Guidewire w/ Trocar Tip, threaded, 0.094 in (2.4 mm) × 9.25 in | AR-8770KT |
| Guidewire w/ Trocar Tip, 0.062 in (1.6 mm) × 7.0 in | AR-8941-7 |

Plates for 6.7 mm/7.0 mm Sets (order separately)

| Product Description | Item Number |
|---|---------------|
| Anterior Plate, 3H, left | AR-8970AL |
| Anterior Plate, 4H, left | AR-8970AL-04 |
| Anterior Plate, 5H, left | AR-8970AL-05 |
| Anterior Plate, 6H, left | AR-8970AL-06 |
| Anterior Plate, 3H, right | AR-8970AR |
| Anterior Plate, 4H, right | AR-8970AR-04 |
| Anterior Plate, 5H, right | AR-8970AR-05 |
| Anterior Plate, 6H, right | AR-8970AR-06 |
| Anterior Plate, short | AR-8970AS-03 |
| Anterior Plate, minimally invasive | AR-8970MA |
| Lateral Tibiotalar Plate, 3H | AR-8970TT |
| Lateral Tibiotalar Plate, 4H | AR-8970TT-04 |
| Lateral Tibiotalar Plate, 5H | AR-8970TT-05 |
| Lateral Tibiotalar Plate, 6H | AR-8970TT-06 |
| Lateral Tibiotalocalcaneal Plate, 3H | AR-8970TTC |
| Lateral Tibiotalocalcaneal Plate, 4H | AR-8970TTC-04 |
| Lateral Tibiotalocalcaneal Plate, 5H | AR-8970TTC-05 |
| Lateral Tibiotalocalcaneal Plate, 6H | AR-8970TTC-06 |
| Posterior Tibiotalocalcaneal Plate, left | AR-8970PL |
| Posterior Tibiotalocalcaneal Plate, right | AR-8970PR |

Low Profile Screws, 4.5 mm/5.5 mm Screws

| Product Description | Item Number |
|---|----------------|
| Low Profile Locking Screws | |
| 4.5 mm × 18 mm-50 mm (2.0 mm increments) | AR-8545L-18-50 |
| 4.5 mm × 55 mm-75 mm (5.0 mm increments) | AR-8545L-55-75 |
| Low Profile Screws | |
| 4.5 mm × 18 mm-50 mm (2.0 mm increments) | AR-8545-18-50 |
| 4.5 mm × 55 mm-100 mm (5.0 mm increments) | AR-8545-55-100 |
| Low Profile Screws, cancellous | |
| 5.5 mm × 20 mm-100 mm (5.0 mm increments) | AR-8555-20-100 |

7.0 XL Compression FT screws

| Product Description | Item Number |
|--|------------------|
| 7.0 XL Compression FT Screws, cannulated, Ti, fully threaded, 40 mm-100 mm (5.0 mm increments) | AR-8770-40H-100H |

Ankle Fusion Plating System, 6.7 mm Set (AR-8970S-67)

| Product Description | Item Number |
|--|----------------|
| Drill Guide, threaded, locking, 4.5 mm, qty. 2 | AR-8970-01 |
| Drill Guide, 3.0 mm/4.5 mm | AR-8970-02 |
| Drill Guide, 4.0 mm/6.7 mm | AR-8967G |
| Drill Guide, 3.0 mm/5.5 mm | AR-8970-05 |
| Depth Measuring Device, long, 4.5 mm/5.5 mm | AR-8970-07L |
| Depth Device, cannulated, for 6.7 mm screws | AR-8967DG |
| Depth Device, large | AR-4167 |
| Drive Shaft, T20 hexalobe, qty. 2 | AR-8970-03 |
| Driver, cannulated, 3.5 mm hex, qty. 2 | AR-8967D |
| Driver, T20 hexalobe, straight | AR-8970-04 |
| Driver, T20 hexalobe, straight, AO, qty.2 | AR-8970-08 |
| Ratcheting Handle, cannulated, large AO handle, QC | AR-8970RH |
| Mini Joint Distractor/Compressor | AR-8970JD |
| Screwdriver Handle, ratcheting | AR-1999 |
| Bone Reduction Forceps, qty. 2 | AR-8943-07 |
| Hohmann Retractor, 9.5 in, 17 mm pointed, qty. 2 | AR-9260-34 |
| Hudson Adapter | AR-1416 |
| Cup Curette, straight shaft, 6.0 mm | AR-8970-11 |
| Cup Curette, curved shaft, 6.0 mm | AR-8970-12 |
| Cobb Elevator, 9.0 mm | AR-8640 |
| Screw Holding Forceps | AR-8941F |
| Countersink, fixed handle, cannulated, 6.7 mm | AR-8967CSF |
| Guidewire Sleeve Insert, 1.6 mm | AR-8970-06 |
| Small Joint Osteotome Angled Up, 0.217 in (5.5 mm) w/ handle | AR-8650-08 |
| Small Joint Osteotome Straight, 0.217 in (5.5 mm) w/ handle | AR-8650-09 |
| Small Joint Osteotome Angled Up, 0.472 in (12 mm) w/ handle | AR-8970-13 |
| Small Joint Osteotome Straight, 0.472 in (12 mm) w/ handle | AR-8970-14 |
| Ankle Fusion Instrument Case | AR-8970C-01 |
| Ankle Fusion Instrument Case, 6.7 mm tray | AR-8970C-67 |
| Ankle Fusion Caddy, 6.7 mm insert | AR-8970C-SC-67 |

References


- Robert KQ, Chandler R, Baratta RV, et. al. The effect of divergent screw placement on the initial strength of plate-to-one fixation. *J Trauma*. 2003;55(6):1139-1144. doi:10.1097/01.TA.0000031103.15337.CA
- Arthrex, Inc. Data on file (APT-05220). Naples, FL; 2021.
- Martin GJ, Boden SD, Titus L, Scarborough NL. New formulations of demineralized bone matrix as a more effective graft alternative in experimental posterolateral lumbar spine arthrodesis. *Spine*. 1999;24(7):637-645. doi:10.1097/00007632-199904010-0000
- Arthrex, Inc. LA1-000006-en-US_A. Naples, FL; 2019.

Disposables for AR-8970S-67 (order separately)

| Product Description | Item Number |
|--|--------------|
| BB-Tak, large | AR-8970-09 |
| BB-Tak, large, threaded | AR-8970-09T |
| Guidewire w/ Trocar Tip, nonthreaded, 0.094 in (2.4 mm) × 8.0 in, qty. 6 | AR-8967K |
| Guidewire w/ Trocar Tip, threaded, 0.094 in (2.4 mm) × 8.0 in, qty. 6 | AR-8967KT |
| Guidewire w/ Trocar Tip, nonthreaded, 0.094 in (2.4 mm) × 12 in, qty. 6 | AR-8967K-12 |
| Guidewire w/ Trocar Tip, threaded, 0.094 in (2.4 mm) × 12 in, qty. 6 | AR-8967KT-12 |
| Guidewire w/ Trocar Tip, 0.062 in (1.6 mm) × 7.0 in, qty. 6 | AR-8941-7 |
| Washer, Ti, 13 mm | AR-8967W |
| Traction Screw, 20 mm | AR-8950JD-2 |
| Drill Bit, cannulated, 4.0 mm | AR-8970-40C |
| Drill Bit, calibrated, long, 3.0 mm | AR-8970-30L |
| Drill Bit, cannulated, long, 3.0 mm | AR-8970-30CL |
| Drill Bit, cannulated, long, 4.5 mm | AR-8970-45CL |
| Drill Bit, long, 4.5 mm | AR-8970-45L |
| Drill Bit, long, 5.5 mm | AR-8970-55L |
| Drill Bit, cannulated, long, 5.5 mm | AR-8970-55CL |

Cannulated Lag Screws

| Product Description | Item Number |
|---|--------------------|
| Low Profile Screws, cannulated, partially threaded, 6.7 mm × 40 mm-100 mm, 18 mm length (5.0 mm increments) | AR-8967-1840-18100 |





This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience, and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level and/or outcomes.

View U.S. patent information at www.arthrex.com/corporate/virtual-patent-marking

arthrex.com

© 2023-07 Arthrex, Inc. All rights reserved. LT1-000197-en-US_B