

Shoulder Suspension System

Quick Setup Guide



1. Attach clamp onto end of table rail
2. Shift patient toward head end on table




3. Insert tower into clamp
4. Pivot forward until secure



5. Disengage pin
6. Rotate boom up
7. Push boom forward



8. Attach weight hanger
9. Push white lever between the handles and lift abduction range arm into desired position



10. Attach arm sleeve connector to sterile disk barrier



11. Attach arm sleeve connector to traction cable pin
12. Rotate blue ring clockwise to tighten (finger tight only; do not overtighten)



13. Attach arm with disposable sleeve
14. Add additional weight



Shoulder Suspension System

Lateral Distraction

15. Attach lateral traction connector to sterile disk barrier
16. Attach connector to lateral boom



17. Extend lateral boom by pulling forward
18. Attach lateral traction sling to connector