
A. Device Description

The Arthrex Shoulder Suspension System (AR-1650S) is an advanced patient positioning system for lateral shoulder arthroscopy. The system permits patient positioning, joint distraction and joint manipulation allowing for arthroscopic treatment of various shoulder conditions. The system is comprised of a Shoulder Suspension Tower, an Arm Sleeve Connector, a Lateral Traction Sling Connector, a Weight Hanger and a set of Weights. A Table Clamp (based on mating table rails) is also included. The system requires the use of the Arm Sleeve Disposables (AR-1651) and a Lateral Traction Sling Disposable (AR-1652) both provided sterile-packaged. A Storage Cart is available to provide safe transport and storage of the entire system. A full listing of components and REF #'s are shown in the tables below.

Item	REF	Description
A	AR-1650-01	Shoulder Suspension Tower
B	AR-1650-02	S3 Arm Sleeve Connector
C	AR-1650-03	S3 Lateral Traction Sling Connector
D	AR-1650-04	Weight Hanger, 10lbs (4.5kg)
E	AR-1650-05	Weight, 2.5lb (1.1kg) – (5 Weights needed per system)
F*	AR-1650-06	Storage Cart, Shoulder Suspension System
G	AR-1650-07	US Clamp, Shoulder Suspension System
H	AR-1651	Disposable, S3 Arm Sleeve
I	AR-1652	Disposable, S3 Lateral Traction Sling

* - Not shown

Market Specific Accessories (one required):

Item	REF	Description
G*	*choose one clamp based on rail std.	
G	AR-1650-08	EU Clamp, Shoulder Suspension System
G	AR-1650-09	UK Clamp, Shoulder Suspension System
G	AR-1650-10	Denyer Clamp, Shoulder Suspension System
G	AR-1650-11	Japan Clamp, Shoulder Suspension System

Product Specifications:

Tower Height (Opened): 62" (157.48cm) from highest point to table top
Cable Travel Range:

- ~26" (66cm) of horizontal travel of cable
- ~40" (102cm) of vertical travel of cable
- 40" (102cm) from center of the clamp to the furthest point of horizontal travel
- Maximum Arm Length: 34" (86cm)
- Maximum Patient Weight: 400lbs (181.44kg)

B. Indications

The Arthrex Shoulder Suspension System is intended to position the patient and apply distal limb suspension and lateral traction during lateral arthroscopic and/or mini-open shoulder procedures.

C. Warnings

- ⚠ Prior to the use, ensure that the system is cleaned and is in proper working condition.
- ⚠ Prior to use, ensure that the connectors (AR-1650-02 and AR-1650-03) have been sterilized. **They are provided NON-STERILE.**
- ⚠ Prior to using this or any other type of medical apparatus with a patient, it is recommended that you read these instructions and familiarize yourself with the product.
- ⚠ Proper understanding of appropriate surgical techniques is required prior to using this product.
- ⚠ Ensure that the patient weight does not exceed the recommended OR table and Shoulder Suspension System specifications.
- ⚠ Ensure that the OR table rails are in good condition. If they are damaged, do not mount the system to the table.
- ⚠ Ensure that the patient's fingers can be folded into a fist as needed to fit into the disposable.
- ⚠ Failure to use this device in accordance with the directions for use may result in device failure, render the device unsuitable for its intended use or compromise the procedure.

D. Precautions

1. Refer to the surgical table manufacturer's user guide for instructions regarding the proper orientation and use of the

Manufactured For:

Arthrex, Inc.
1370 Creekside Blvd
Naples, FL 34108, USA
1-(800) 934-4404

Arthrex GmbH
Erwin-Hielscher-Strasse 9
81249 München, Germany
+49 89 909005-0

Allen Medical Systems
100 Discovery Way
Acton, MA 01720
(978) 263-5401

Hill Rom SAS
B.P. 14 Z.I. du Talhouet
56330 Pluvigner, France
Tel. +33 (0) 2 97 50 92 12
Fax. +33 (0) 2 97 50 92 12

table. Always refer to the surgical table manufacturer’s weight limits.

2. Ensure that the OR table is locked prior to mounting the Shoulder Suspension System.
3. Do not use if product shows visible damage.
4. Ensure device is securely attached and locked to OR Table accessory rail before transferring patient.
5. Different OR tables may have different table rail dimensions depending on country of manufacture/use. Ensure that the correct Table Clamps designed for use with the specific OR table rail are used.
6. Pinch points and crushing points are indicated on the device via labeling shown in the following section.

E. Labeling

	Single Use only: this symbol means that the product cannot be reused.
	Crush force: this symbol indicates that a pinching force can occur in this area. Use Caution.
	Pinch force: this symbol indicates that a pinch point is present. Use Caution.
	Do Not Immerse: This symbol means that the product should not be submerged in water or cleaning solutions.
	Home: this symbol indicates the lowest position for the Pulley Carriage Bars.
	One: this symbol indicates the middle position for the Pulley Carriage Bars
	Two: this symbol indicates the highest position for the Pulley Carriage Bars
	Maximum Traction Force: this symbol indicates the highest amount of force that can be applied


Push Prohibited: Pushing this item without using the designated hold positions can result in hazardous situation

F. Equipment Setup

- i** Note: Ensure that the surgical table is in its final operation location with respect to the room.
- i** Note: Remove X-ray tops if present from the surgical table.
- i** Note: Ensure that the longest part of the surgical table base is located under the Suspension System attachment location for greatest table stability.
- i** Note: Ensure that the OR table wheel locks are activated before placing the patient on the system.

Clamps (AR-1650-07, AR-1650-08, AR-1650-09, AR-1650-10, or AR-1650-11)

Select the appropriate Clamp (refer to Item G) based on the OR table and table rail dimensions that the system will be mounted to. Each rail clamp is designed for specific rail dimensions. Use the following diagram and table to determine the correct clamp that should be used with your table and rail type.

Rail Type	Height	Width
American	1 1/8"	3/8"
European	25mm	10mm
United Kingdom	1 ¼"	¼"
Denyer	1 ½"	¼"
Japan	32mm	9mm

Manufactured For:

Arthrex, Inc.
1370 Creekside Blvd
Naples, FL 34108, USA
1-(800) 934-4404

Arthrex GmbH
Erwin-Hielscher-Strasse 9
81249 München, Germany
+49 89 909005-0

Allen Medical Systems
100 Discovery Way
Acton, MA 01720
(978) 263-5401

EC REP

Hill Rom SAS
B.P. 14 Z.I. du Talhouet
56330 Pluvigner, France
Tel. +33 (0) 2 97 50 92 12
Fax. +33 (0) 2 97 50 92 12

Figure 1

Mount Clamp to the rail by sliding it on the distal end of the rail. Typically, the Clamp is placed on the contralateral side of the patient but it is not required. It is recommended for most tables and patients that the clamp be pulled back to engage the rail gravity lock. Secure the Table Clamp to the rails by tightening the two wing knobs on the Clamp body (Figure 1).

Lift the lever, move it approximately to center and push the lever down (Figure 2). Adjusting the lever forward or backward will adjust the angle of the Tower relative to the OR table. It is recommended that the starting orientation is perpendicular to the rail in order to mount the Tower parallel with the patient but it is not required.

Figure 2

Suspension Tower

Remove the Suspension Tower from the Storage Cart or designated storage area by holding it by the blue carrying handles shown in Figure 3.

Figure 3

Place the hook end of the Tower (AR-1650-01) into the slot in the stem of the clamp as shown (Figure 3). Lift and push forward on the Tower to pivot it in the Clamp (Figure 4). The Tower will rest at about 15 degrees toward the head of the patient.

Manufactured For:

Arthrex, Inc.
1370 Creekside Blvd
Naples, FL 34108, USA
1-(800) 934-4404

Arthrex GmbH
Erwin-Hielscher-Strasse 9
81249 München, Germany
+49 89 909005-0

Allen Medical Systems
100 Discovery Way
Acton, MA 01720
(978) 263-5401

Hill Rom SAS
B.P. 14 Z.I. du Talhouet
56330 Pluvigner, France
Tel. +33 (0) 2 97 50 92 12
Fax. +33 (0) 2 97 50 92 12

Figure 4

When in place, the blue latch on the Tower hook will snap into place, preventing the Tower from being inadvertently tipped back. If the surgical table pad is askew, it may touch the Clamp and interfere with the latch. Move the table pad or push the latch into place.

 Caution: Ensure that latch has been fully engaged before proceeding

Figure 5

To unfold the Tower, pull the black knob adjacent to the Tower hook (Figure 6) to unlock the Tower. Lift up on the handle to pivot the boom of the Tower (Figure 7).

Figure 6

Push the Tower boom forward until it clicks into place. (Figure 7)
The kickstand should also be dropped into place in the orientation shown.

Manufactured For:

Arthrex, Inc.
1370 Creekside Blvd
Naples, FL 34108, USA
1-(800) 934-4404

Arthrex GmbH
Erwin-Hielscher-Strasse 9
81249 München, Germany
+49 89 909005-0

Allen Medical Systems
100 Discovery Way
Acton, MA 01720
(978) 263-5401

Hill Rom SAS
B.P. 14 Z.I. du Talhouet
56330 Pluvigner, France
Tel. +33 (0) 2 97 50 92 12
Fax. +33 (0) 2 97 50 92 12

Figure 7

 Caution: Upon unfolding Tower, inspect cable for any fraying, kinking or damage

Weight Hanger

Weight Hanger (AR-1650-04) attaches via its hook to the triangular receiver which hangs from the back of the Tower (Figure 8). The Weight Hanger weighs 10 lbs (4.5kg).

Figure 8

Weights

Each Weight (AR-1650-05) weighs 2.5lbs (1.1kg). They are stacked onto the Weight Hanger by aligning the slot of the Weight with the post of the Weight Hanger and sliding down over the cone of the Weight Hanger (Figure 9).

Figure 9

G. Sterile Setup

Positioning the Patient

Ensure the patient is laying in the lateral position with the operative shoulder facing upward.

Manufactured For:

Arthrex, Inc.
1370 Creekside Blvd
Naples, FL 34108, USA
1-(800) 934-4404

Arthrex GmbH
Erwin-Hielscher-Strasse 9
81249 München, Germany
+49 89 909005-0

Allen Medical Systems
100 Discovery Way
Acton, MA 01720
(978) 263-5401

Hill Rom SAS
B.P. 14 Z.I. du Talhouet
56330 Pluvigner, France
Tel. +33 (0) 2 97 50 92 12
Fax. +33 (0) 2 97 50 92 12

Positioning the Tower

Position the Tower so that it is oriented over the patient's operative arm. The angle of the Tower can be adjusted by releasing the trigger and rotating it either forward or backward. This should be done by a circulating nurse or other qualified non sterile personnel.

Adjust the height of the cable as needed. This can be done by pressing the trigger and moving the pulley carriage to any of the three main ranges (Figure 10). Note that the trigger will lock at each of the three ranges and require activation to move further.

Figure 10

Arm Sleeve Connector

The Arm Sleeve Connector is comprised of three different sections. The top section in blue can be rotated to loosen or tighten the grip of the hook at the top (Figure 11). The bottom section in black is where the Arm Sleeve will click into place. The grey bell section can be rotated in order to lock the black clip section in place to restrict rotation of the arm. Finally there is a plastic circle divider separating the sterile and non-sterile fields. The plastic divider can be found in the Arm Sleeve Disposable package.

Figure 11

The disposable Arm Sleeve Connector should be mounted in the pin at the center of the traction cable. The hook used to attach the Arm Sleeve Connector can be loosened (by twisting the blue top section clockwise) and tightened (by twisting the blue top section counter clockwise) once mounted to engage the pin.

Loosen the lower section of the Connector (by twisting the grey bell section clockwise) in order to allow for the patients arm to freely rotate during initial setup. Note the location of the plastic separating the sterile field.

Manufactured For:

Arthrex, Inc.
1370 Creekside Blvd
Naples, FL 34108, USA
1-(800) 934-4404

Arthrex GmbH
Erwin-Hielscher-Strasse 9
81249 München, Germany
+49 89 909005-0

Allen Medical Systems
100 Discovery Way
Acton, MA 01720
(978) 263-5401

Hill Rom SAS
B.P. 14 Z.I. du Talhouet
56330 Pluvigner, France
Tel. +33 (0) 2 97 50 92 12
Fax. +33 (0) 2 97 50 92 12

Arm Sleeve Disposable (AR-1651)

Once the patient is on the table and draped, locate the Arm Sleeve disposable and wrap it around the patient's arm in the manner prescribed and shown in the instructional artwork on the disposable. (Figure 12)


Figure 12

⚠ Use only the approved disposable, AR-1651 Arm Sleeve, Non-conformance may result in patient not being secured properly

First unwrap the disposable so that it can lay flat. Place the patients hand in the sleeve so that it rests below the white clip and on top of the flat section with numbered edges 1-4 on the reverse side. Wrap sections 1 & 2 and then follow with 3 & 4 (Figure 13). Flip the second flat section with numbered edges 5-7 over the patients arm so that it is covered and the white clip is facing out away from the patients knuckles. Finally secure straps 5, 6 and 7 (Figure 14).

Figure 13

Use the numbered tabs to show the order in which segments of the disposable should be wrapped on the patient's operative arm.

Figure 14

Once the arm sleeve is wrapped around the operative arm, pull the Arm Sleeve Connector to the arm sleeve by grabbing it by the grey bell section of the connector and click the arm sleeve in place (Figure 15). Ensure that the two white clips are exposed so that the patient is fully locked in. Apply traction to the patient's arm as needed by adding or removing weight from the weight hanger. This should be performed by a circulating nurse or other qualified non sterile personnel.

Figure 15

Adjust rotation of the arm as needed. The arm can be held in the desired location by tightening (counterclockwise rotation) the grey bell on the bottom of the Arm Sleeve Connector.

To adjust the angle at which the arm is held, grab the connector by the grey bell section and slide to the desired location (Figure 16).

Manufactured For:

Arthrex, Inc.
 1370 Creekside Blvd
 Naples, FL 34108, USA
 1-(800) 934-4404

Arthrex GmbH
 Erwin-Hielscher-Strasse 9
 81249 München, Germany
 +49 89 909005-0

Allen Medical Systems
 100 Discovery Way
 Acton, MA 01720
 (978) 263-5401

EC REP
 Hill Rom SAS
 B.P. 14 Z.I. du Talhouet
 56330 Pluvigner, France
 Tel. +33 (0) 2 97 50 92 12
 Fax. +33 (0) 2 97 50 92 12

Figure 16

Figure 18

Mount the Lateral Traction Sling by attaching the loop end of the Lateral Traction Sling on the hook on the bottom of the Lateral Traction Connector (Figure 19).

Lateral Sling Connector

Attach the Lateral Sling Connector by hooking it to the end of the Tower furthest from the clamp (Figure 17).

Figure 17

Once the hook has grabbed the pin, pull the lateral suspension arm into place using the Lateral Sling Connector. The Lateral Sling's position can be manipulated by adjusting the position of the Connector.

Lateral Traction Sling Disposable (AR-1652)

Put the patient's arm in the sling so that the patient's arm is being held up by the Lateral Traction Sling (Figure 18).

Figure 19

Apply lateral traction by reducing the length of the Lateral Traction Sling. The length of the sling can be adjusted by either pulling down on the loop to increase the traction applied or reduced by pressing on the white clip. Use the tensiometer on the Lateral Traction Sling connector to determine the amount of traction needed to properly position the patient (Figure 20).

Caution: Do not exceed 30lbs (13.6kg) traction

Manufactured For:

Arthrex, Inc.
 1370 Creekside Blvd
 Naples, FL 34108, USA
 1-(800) 934-4404

Arthrex GmbH
 Erwin-Hielscher-Strasse 9
 81249 München, Germany
 +49 89 909005-0

Allen Medical Systems
 100 Discovery Way
 Acton, MA 01720
 (978) 263-5401

Hill Rom SAS
 B.P. 14 Z.I. du Talhouet
 56330 Pluvigner, France
 Tel. +33 (0) 2 97 50 92 12
 Fax. +33 (0) 2 97 50 92 12

Figure 20

The Lateral Suspension Arm can be moved anteriorly or posteriorly in regards to the patient. This may be helpful to apply an opposite traction force to the shoulder joint. The Lateral Sling's tension must be released prior to moving the Lateral Suspension Arm. (Figure 21)

Figure 21

H. Removal and Storage

Removal of Disposable Components

Lateral Traction Sling

Remove the Lateral Traction Sling from the hook at the end of the Lateral Sling Connector. Dispose of sling after use.

Arm Sleeve

Remove the Arm Sleeve from the Arm Connector by pressing down on the white buttons to unclip the sleeve from the connector (Figure 22). Dispose of arm sleeve after use.

Ensure that you are supporting the patient's arm before releasing the sleeve from the connector.

Figure 22

Manufactured For:

Arthrex, Inc.
1370 Creekside Blvd
Naples, FL 34108, USA
1-(800) 934-4404

Arthrex GmbH
Erwin-Hielscher-Strasse 9
81249 München, Germany
+49 89 909005-0

Allen Medical Systems
100 Discovery Way
Acton, MA 01720
(978) 263-5401

Hill Rom SAS
B.P. 14 Z.I. du Talhouet
56330 Pluvigner, France
Tel. +33 (0) 2 97 50 92 12
Fax. +33 (0) 2 97 50 92 12

Removal of System Components

Lateral Sling Connector

Remove the Lateral Sling Connector by pressing up on the silver trigger to release the lock. Then lift the connector off of the lateral extension rod (Figure 23).

Figure 23

Arm Sleeve Connector

Remove the Arm Sleeve Connector by loosening the top blue section of the connector (clockwise). Once the hook is loose remove the connector.

Suspension Tower

To collapse the Tower first ensure that the pulley carriage is in the Home position (Figure 24).

Figure 24

Once the pulley carriage is in the bottom Home position pull down on the handle at the end of the Tower to release the lock. Allow the top of the Tower to slide back towards the operator until the kickstand disengages (Figure 25).

Figure 25

Once the kickstand has disengaged press the upper end of the Tower where the handle is located into its Home section on the lower section of the Tower. There should be a click when the Tower is completely folded and locked in place (Figure 26).


Figure 26

Remove the Tower by pressing up on the blue trigger at the bottom of the Tower. This allows the Tower to be tilted back and removed. Place the Tower in the Cart (AR-1650-06).

Manufactured For:

Arthrex, Inc.
1370 Creekside Blvd
Naples, FL 34108, USA
1-(800) 934-4404

Arthrex GmbH
Erwin-Hielscher-Strasse 9
81249 München, Germany
+49 89 909005-0

Allen Medical Systems
100 Discovery Way
Acton, MA 01720
(978) 263-5401

Hill Rom SAS
B.P. 14 Z.I. du Talhouet
56330 Pluvigner, France
Tel. +33 (0) 2 97 50 92 12
Fax. +33 (0) 2 97 50 92 12

Table Clamps

Loosen the knobs on the lateral aspect of the Clamp body to unlock the Clamp from the table rail. Slide the Table Clamp off the rail, disengaging the table rail gravity locks if appropriate. Clean as directed in the following section and store in a safe location, preferably the designated location on the Storage Cart.

Storage Cart (AR-1650-06)

Careful storage of the Arthrex Shoulder Suspension System can help prolong its life. It is recommended that the system be stored in the Storage Cart (AR-1650-06). The Cart features locations for each component to secure and protect them and make transportation easier.

The Storage Cart can be loaded as shown (Figure 27). Lower the Tower into its place in the Cart by placing the hook end of the Tower in the slot between the cart body and the pin. To store the Table Clamp, slide it over the end of the Cart latch and lower the latch into place. The Cart will hold up to six Weights, the Weight Hanger, the Tower and the Clamp. The Cart allows for simple movement and storage of the system. The Cart should be stored in a location with an even surface.

Figure 27

⚠ Caution: Do not leave Cart on a inclined surface

⚠ Warning: Applying force to the Tower when it is loaded in the Cart can result in tipping the cart

I. Cleaning and Maintenance

After each use, clean the device. Follow the instructions for cleaning of the applicable parts below

⚠ Caution: Do not clean with bleach. Failure to properly clean product can result in reduced performance.

- Tower, Weights, Cart and Clamps:
 - Clean and disinfect all surfaces using a quaternary ammonium disinfecting/cleaning solution following the manufacturer’s recommendation for achieving low-level disinfection.
 - Read the cleaning product’s directions and follow the instructions on the label.
 - Use caution in areas where fluid migration may occur.
 - Wipe device with a clean, dry cloth.
 - Make certain the product is dry prior to reinstalling to avoid damage.
- Sterile Connections:
 - Refer to Sterilization and Cleaning instructions (D-720771)
- Disposables
 - Dispose of materials appropriately according to facility requirements. **DO NOT USE MULTIPLE TIMES.**

Should any need for maintenance or repair occur please contact customer service, contact information below.

J. Customer Service Contact Information

	North America	Latin America	Europe and other Countries
☎	1-866-267-9138 1-239-643-5553, x5603	001 954-447-6815	0049 8131 5957 0
📠	1-239-591-6943	001 954-447-6814	0049 8131 5957 631

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