

Virtual Implant Positioning™

VIP™ Reamer

Depth of Ream Instrumentation with Positive Stop
for Confident Precise Reaming



- Reusable patient-specific instrumentation designed to reduce cost and lead time
- Confidently transfer your planned depth of ream with positive-stop technology
- Planned and executed reaming in 1 mm increments for VaultLock and MGS glenoid systems
- Pilot reamer creates a flat bottom hole followed by a secondary reamer



Learn more about the
VIP reamer



© 2023-08 Arthrex, Inc. All rights reserved.
FL1-000694-en-US_A