









# Reference Chart

## MIS Burrs

Below is a quick reference guide for the Arthrex MIS burrs.

Each MIS burr is designed for a specific procedural application and it is imperative that the proper burr is used. Please refer to the chart below to select the appropriate burr for each procedure.

	Burr	Product Description	Indication	Item Number
Hallux Valgus Correction		2.0 mm × 13 mm, straight	Akin Osteotomy	AR-300-B001
		2.0 mm × 19.5 mm, straight	1st Metatarsal Osteotomy (Chevron and Transverse)	AR-300-B201
Calcaneal Osteotomy		3.1 mm × 20 mm, straight	Calcaneal Osteotomy	AR-300-B202
Cheilectomy		4.3 mm × 13 mm, conical – wedge	Cheilectomy or Bone Resection	AR-300-B101
		2.9 mm × 13 mm, straight – wedge	Cheilectomy or Bone Resection	AR-300-B102
		5.0 mm × 15 mm, oval	Large Bone Resection	AR-300-B103
Lesser Toe Osteotomy		2.0 mm × 8 mm, straight	Hammertoe Correction	AR-300-B002
		2.2 mm × 12 mm, straight	Tailor's Bunion	AR-300-B003