


# Virtual Implant Positioning<sup>™</sup> Reamer

Reusable Transfer Technology with Positive Stop for  
Precise Reaming



- Reusable patient-specific instrumentation designed to reduce cost and lead time
- Planned and executed reaming in 1 mm increments for Univers VaultLock<sup>®</sup> and MGS glenoid systems
- Confidently transfer your planned depth of ream with positive-stop technology
- Pilot reamer creates a flat bottom hole followed by a secondary reamer



Learn more about the  
VIP<sup>™</sup> reamer



© 2024-10 Arthrex, Inc. All rights reserved.  
FL1-000694-en-US\_B