


SwiftStitch™ Suture Passer

Product Highlights



Features and Benefits

Low Profile

- The 1.5 mm outer diameter creates a minimal hole through the labrum

Versatile

- Accommodates suture sizes ranging from #2-0 FiberWire® suture to SutureTape


Multifunctional

- Can be used for simple or mattress suture configurations, and both knotted and knotless labral repair techniques

Precise

- Button control of the suture allows surgeons to release it from the jaws in desired location

Retrieve and Release



Ordering Information

Product Description	Item Number
SwiftStitch Suture Passer, Crescent Shape	AR-4068HL

Products may not be available in all markets because product availability is subject to the regulatory approvals and medical practices in individual markets. Please contact Arthrex if you have questions about the availability of products in your area.

View U.S. patent information at www.arthrex.com/corporate/virtual-patent-marking

© 2018 Arthrex, Inc. All rights reserved. FL1-00107-EN_A

