

The Foundation of Research-Based Shoulder Arthroplasty Innovation


Advance with us...




Eclipse™

Total Shoulder System: The Cage-Screw Difference

- 14+ years of long-term survivorship data¹
- Proprietary cage-screw fixation provides dynamic compression at the bone-prosthesis interface and reduces medial calcar osteolysis²
- Reproducible SpeedScap™ subscapularis repair system
- Extensive glenoid augment options for treating simple and complex pathologies
- Single instrument tray allows streamlined three-step implantation
- Revision-friendly design






REVERSE

Univers Revers™

Total Shoulder System

- 135° neck-shaft angle optimizes impingement-free range of motion^{3,4}
- Inlay humeral design facilitates ideal soft-tissue tensioning⁵
- Fracture repair features include press-fit stem geometry, 135° inclination angle, and soft-tissue fixation points for improved functional outcomes and reliable tuberosity healing⁶
- Designed-to-fit Modular Glenoid System includes screw or post fixation, augments, and multiple lateralization and glenosphere options






VIP

The **VIP**TM Experience





For Patients

- MyVIP delivers timely, automated text messages to guide patients through their surgical journey
- Clinic tools create valuable opportunities for educating patients and easing anxiety-ridden waiting time
- Curated web and social media content help with patient outreach





For Surgeons

- Cross-platform capability across all devices
- Improves surgical decision-making and implant positioning⁷⁻¹⁰
- 24- to 48-hour turnaround after plan approval
- Precisely execute your planned depth of ream with reusable patient-specific instrumentation
- Reusable patient-specific instrumentation eliminates excessive costs and lead time¹¹
- Optimize your humeral component selection to achieve maximum fixation in any anatomy

Humeral Planning Now Available



Global Leaders in Surgical Skills Education

References

1. Habermeyer P, et al. *JSES Open Access*. 2019;3(4):234. doi:10.1016/j.jses.2019.10.013
2. Alikhah A, et al. *BMC Musculoskeletal Disord*. 2020;21(1):295. doi:10.1186/s12891-020-03277-3
3. Walch G, *J Shoulder Elbow Surg*. 2019;28(6):1214-1216. doi:10.1016/j.jse.2019.02.018
4. Gobezie R, et al. *J Shoulder Elbow Surg*. 2019;28(5):813-818. doi:10.1016/j.jse.2018.11.064
5. Haidamous G, et al. *J Shoulder Elbow Surg*. 2020;29(12):2556-2563. doi:10.1016/j.jse.2020.03.036
6. O'Sullivan, et al. *J Shoulder Elbow Surg*. 2020;29(9):1938-1949. doi:10.1016/j.jse.2020.03.032
7. Werner BS, et al. *J Shoulder Elbow Surg*. 2017;26(8):1477-1483. doi:10.1016/j.jse.2017.01.006
8. Iannotti JP, et al. *J Bone Joint Surg Am*. 2019;101(5):446-457. doi:10.2106/JBJS.17.01614
9. Heylen S, et al. *J Shoulder Elbow Surg*. 2016;25(2):186-192. doi:10.1016/j.jse.2015.07.024
10. Scalise JJ, et al. *J Bone Joint Surg Am*. 2008;90(11):2438-2445. doi:10.2106/JBJS.G.01341
11. Cancienne JM, et al. *Shoulder Elbow*. 2021;13(4):380-387. doi:10.1177/178573220908865



Scan for more information
about Arthrex Arthroplasty
innovations!

