


ARTHREX

PIVOTPOST



The PivotPost is a lateral post that moves in three planes to easily place in the proper position to the patient, and pivot out of the way when not in use.

PIVOT



IN/OUT




UP/DOWN

FEATURES AND BENEFITS:

- Prevent patient's leg from moving over post when valgus stress is applied
- Three planes of motion (rotation, in/out, up/down)
- Small size allows surgeons to operate around easily
- Minimal weight and simple set-up for OR staff
- Disposable pad that is easy to put on and take off
- Ergonomic, oversized buttons are easily located under the drape



The unique over-the-top sweeping design prevents the patient's leg from moving over the post when valgus stress or a figure four is applied.



CLARK RAIL ADAPTER

- Custom for the PivotPost
- Stainless steel
- Finger tight holds firmly
- The Pull Knob to raise and lower is round
- The Tightening Knob is scalloped



FOAM PAD

- Disposable (nonsterile)
- Five in a box
- No specific top or bottom
- Double-sided Velcro®



ORDERING INFORMATION

PivotPost	AR-1560
Knee Post only	AR-1560-02
Clark Rail Adapter only	AR-1560-03
Foam Pads (box of five)	AR-1560-01

www.arthrex.com

... up-to-date technology
just a click away

U.S. PATENT PENDING

©2014, Arthrex Inc. All rights reserved. LS1-0165-EN_A